



CUBIC-HV™1 3D staining protocol (version 2020.07)

[for a whole mouse brain]

[Technical note - c-Fos staining](#)

CUBIC-HV™1 3D nuclear staining kit (#C3698)

CUBIC-HV™1 3D immunostaining kit (#C3699)



CUBICStars
Make the world clear

© 2020 CUBICStars Co., All rights reserved

MATERIALS

Reagents in kit

CUBIC-HV™1 3D nuclear staining buffer (1x) (Store at room temperature)

CUBIC-HV™1 3D immunostaining buffer (2x) (Store at 4°C; use within 6 months)

CUBIC-HV™1 3D immunostaining wash buffer (1x) (Store at room temperature; cooled to 4°C when in use)

CUBIC-HV™1 3D immunostaining additive (10x) (included in CUBIC-HV™1 3D immunostaining kit; Store at room temperature, protected from light)

Other reagents

PBS (Tablet): TaKaRa #T9181

HEPES: Tokyo Chemical Industry #H0396

Sodium azide (NaN₃): nacalai tesque #31208-82

Paraformaldehyde: nacalai tesque #02890-45

Formalin solution: Tokyo Chemical Industry #F0622

Heparin: FUJIFILM Wako #081-00136

CUBIC-L: Tokyo Chemical Industry #T3740

CUBIC-R+: Tokyo Chemical Industry #T3741

Nuclear staining reagents

DAPI (included in CUBIC-HV™1 3D nuclear staining kit)

BOBO™-1 Iodide: ThermoFisher Scientific B3582

SYTOX™ Green Nucleic Acid Stain: ThermoFisher Scientific S7020

Propidium Iodide (PI) (included in CUBIC-HV™1 3D nuclear staining kit)

RedDot™2 Far-Red Nuclear Stain: Biotium #40061

Antibodies

Primary antibody

Anti-c-Fos: CST #2250S

[NOTE] To use this antibody, a custom product with a concentration of 200 ug/mL or higher is required.

Secondary antibody

FabuLight Fc specific Fab fragment (Jackson Immunolab)

<https://www.jacksonimmuno.com/catalog/31#target:15>

Fc-Fab-A594: Jackson Immuno Research #111-587-008

Fc-Fab-Cy3: Jackson Immuno Research #111-167-008

Fc-Fab-A647: Jackson Immuno Research #111-607-008

[NOTE] Alexa Fluor™ 488 is not compatible with CUBIC-R+. Cy3, Alexa Fluor™ 594 and 647 have been validated.

Containers to be used

Protein LoBind 500 μ L tube: Eppendorf #022431064

5 mL tube (included in CUBIC-HVTM1 3D nuclear staining kit)

15 mL standing tube (included in CUBIC-HVTM1 3D immunostaining kit)

30 mL tube: SARSTEDT #60.544

50 mL tube: Falcon #352070

Preparation of stock reagents

3D nuclear staining wash buffer (10 mM HEPES, pH7.5)

1M HEPES (pH 7.5)	5 mL
-------------------	------

Distilled water	495 mL
-----------------	--------

=====

Total 500 mL supplied with 0.05% NaN₃

Process outline

Perfusion fixation and dissection of the mouse brain

↓

Post-fix

↓ ~24 h

Wash (PBS)

↓ 2 h x3

Delipidation (CUBIC-L)

↓ 3~7 days

Wash (PBS)

↓ 2 h x3

Nuclear staining (HV™1 3D nuclear staining buffer)

↓ 3~5 days

Wash (3D nuclear staining wash buffer)

↓ 2 h x3

Primary antibody + secondary antibody reaction / Replacement with HV™1 3D immunostaining buffer

↓ 1.5 h

Immunostaining (HV™1 3D immunostaining buffer)

↓ 2-3 weeks

Reaction at 4°C

↓ 1 day

Wash (HV™1 3D immunostaining wash buffer)

↓ 2 h x3

Post-fix

↓ 1 day

Wash (PBS)

↓ 2 h

RI matching (CUBIC-R+)

↓ ~3 days

Microscopic observation

Process details

[NOTE] Reagent volume and reaction time are indicated for staining a single whole mouse brain. It is required to adjust the volume and reaction time according to sample size.

[NOTE] All shaking steps except for enzyme reaction, immunostaining, and 4°C reaction are performed with the tube in the horizontal position.

[NOTE] Brain samples after delipidation are easily damaged and should be handled with a metal spoon.

1) Collection of mouse brain*

1. Anesthetize the mouse with an overdose of pentobarbital sodium salt (nacalai tesque #02095-04) in PBS or saline.
2. Transcardially perfuse with 10 mL (4 mL/min) of cold heparin-PBS (+ 10 U/mL Heparin).
3. Transcardially perfuse 20 mL (6 mL/min) of cold 4% (w/v) paraformaldehyde.
4. Dissect the brain from the skull.
5. Post-fix the dissected brain in 4% (w/v) PFA in PBS (~10 mL/whole brain) for overnight (8-24 h) at 4°C with gentle shaking (40-50 rpm/min).
6. Wash the sample in PBS (+0.05% NaN₃) for 3 h x 3 times at room temperature with gentle shaking (40-50 rpm/min).

*Refer to Susaki et al. Nature Protocol 10:1709–1727 (2015) for details.

[NOTE] Prolonged fixation or storage may result in inadequate detection of c-Fos. It is recommended to proceed to step 2) immediately after PBS wash.

2) Delipidation

1. Immerse the fixed sample in 10 mL of 0.5x CUBIC-L (1:1 dilution with distilled water) in the 30 mL tube and incubate it for overnight at room temperature with gentle shaking (40-50 rpm/min).
2. Replace to 10-15 mL of 1x CUBIC-L in the 30 mL tube and delipidate for 3-5 days* at 37°C with gentle shaking (40-50 rpm/min).

*If the CUBIC-L treatment is more than 3 days, replace with a new CUBIC-L every 2 to 3 days.

3. Wash the sample with 20 mL of PBS (+0.05% NaN₃) for 2 h x 3 times (or 2 h x1, overnight x1, 2 h x1) at 37°C with gentle shaking (40-50 rpm/min).

[NOTE] The tubes should be washed or replaced each time to intensively remove Triton X-100.

[NOTE] The delipidated sample can be stored in PBS/NaN₃ at 4°C.

3) 3D nuclear staining

1. Dilute either of nuclear stains in 3-4 mL of 1x HVTM1 3D nuclear staining buffer in the 5 mL tube.

DAPI: 1/200

BOBO-1: 1/400

SYTOX-G: 1/2500

PI: 1/100

RedDot2: 1/250

5) 3D Immunostaining

1. Prepare antibody staining solution (500 μ L per whole brain) as follows:

2x HV TM 1 3D immunostaining buffer	250 μ L (final 1x)
10x HV TM 1 additive	5 μ L (final 0.1x)
Distilled water	245-(X+Y) μ L

=====

Add the entire volume to the 500 μ L tube containing an X+Y μ L antibody mixture.

2. Transfer the mixed antibody staining solution into the 15 mL standing tube and immerse the buffer-exchanged sample. Incubate it protected from light for 3 weeks at 32°C with gentle shaking (40-50 rpm/min).

[NOTE] Shake the tube in an upright position as during the staining.

[NOTE] To avoid damage, put the brain so that the dorsal side comes to the bottom of the tube.

[NOTE] Close the lid tightly and wrap the parafilm to prevent evaporation.

3. To stabilize the 2nd Fab signals, gently shake the staining tube protected from light with a shaker (40-50 rpm/min) for a further 24 h at 4°C.

[NOTE] Shake the tube in an upright position as during the staining.

4. Wash the sample in 15 mL of pre-cooled 1x HVTM1 3D immunostaining wash buffer in the 30 mL tube protected from light for 30 min x 2 times at 4°C with gentle shaking (40-50 rpm/min).

[NOTE] To stabilize the binding of Fab antibodies and reduce the non-specific signal, make the reagents cooled to 4°C prior to washing and post-staining fixation operations.

6) Post-fixation

1. Prepare the fixative solution by diluting the saturated formalin (FA) product to the final 1% in 1x HVTM1 3D immunostaining wash buffer and cool it to 4°C.

[NOTE] The saturated formalin solution contains 35 to 38% formalin. For example, when you use a 37% formalin solution, dilute it with 1x HVTM1 3D immunostaining wash buffer at a ratio of 1:36.

2. Immerse the sample in 8 mL of 1% FA in the 15 mL standing tube protected from light for 24 h at 4°C with gentle shaking (40-50 rpm/min).

3. To accelerate the fixation reaction, further incubate the sample in 1% FA protected from light for 1 h at 37°C with gentle shaking (40-50 rpm/min).

4. Wash the sample in 15 mL of PBS in a 30 mL tube protected from light for 2 h at 25°C with gentle shaking (40-50 rpm/min).

7) RI matching

1. Immerse the sample in 15 mL of 0.5x CUBIC-R+ (1:1 diluted with water) in the 30 mL tube protected from light for 24 h at 25°C with gentle shaking (40-50 rpm/min).

2. Exchange with 15 mL (or 30 mL if gel embedding is performed) of 1x CUBIC-R+ for 2 days at 25°C with gentle shaking (40-50 rpm/min).

3. Use for microscopic observation (embed the sample in gel if necessary).

REFERENCES

1. Susaki et al. Versatile whole-organ/body staining and imaging based on electrolyte-gel properties of biological tissue. *Nature Communications* (2020) 11: 1982. DOI: 10.1038/s41467-020-15906-5
2. Matsumoto et al. Advanced CUBIC tissue clearing for whole-organ cell profiling. *Nature Protocols* (2019) 14: 3506–3537. DOI: 10.1038/s41596-019-0240-9
3. Susaki et al. Advanced CUBIC protocols for whole-brain and whole-body clearing and imaging. *Nature Protocols* (2015) 10: 1709–1727. DOI: 10.1038/nprot.2015.085

ORDERING INFORMATION

CUBIC-HV™1 3D nuclear staining kit (Tokyo Chemical Industry #C3698)
CUBIC-HV™1 3D immunostaining kit (Tokyo Chemical Industry #C3699)

RELATED PRODUCTS

CUBIC-L (Tokyo Chemical Industry #T3740)
CUBIC-R+ (Tokyo Chemical Industry #T3741)
Formalin solution (Tokyo Chemical Industry #F0622)
Mounting Solution (RI 1.520) [for CUBIC-R+] (Tokyo Chemical Industry #M3294)

CONTACT US



Company website: <https://www.cubicstars.com>

Product website: <https://www.cubicstars.com/cubic-hv/index.html>

CUBIC-HV™ support by email: cubic-hv@cubicstars.com

CUBIC-HV™ support via twitter: [@cubic_hv](https://twitter.com/cubic_hv)